


DATE	SYM	REVISION RECORD	AUTH.	DR.	CK.
08.00	A	Added pin 2 to IL4F	MJ	CB	



	P1 AG-206	P2 AG-206	P3 IL4F	SIGNAL
PIN 1	PIN 1	PIN 1	PIN 1,2	POWER -
PIN 2	PIN 2	PIN 2	---	V0 OXY CURRENT
PIN 3	---	---	---	---
PIN 4	PIN 4	PIN 4	---	V1 OXY TEMP
PIN 5	---	---	PIN 3	V2 pH
PIN 6	PIN 6	PIN 6	PIN 4	POWER +

PART NUMBER	A	B	C
171432	25 IN	28 IN	25 IN

TOLERANCES	SEA-BIRD ELECTRONICS, INC			
FRACTIONAL	P/N SEE TABLE	SCALE NTS	DRAWN BY MJ	
DECIMAL	TITLE Y-CABLE: SBE23 & AMT pH			
ANGULAR	DATE 08.02.00	DRAWING NUMBER	32435	REV A